

MACH®

Precise Working Dies in Just Minutes

- ▶ Flows beautifully—better than stone—to create a smoother die and to capture minute detail.
- ▶ No hand-mixing—inject it directly from the automix cartridge into impression.
- ▶ Exceptionally hard set—90+ Durometer.
- ▶ Slightly flexible.
- ▶ Set time: 2 minutes



SKU	DESCRIPTION
S433S	Mach® Die-Silicone Includes (2) 50ml automix cartridges (100ml ttl.), (10) Yellow Mixing tips, and (10) Dark Grey Intraoral tips

Accessories & Refills

S302	50 ml Mixing Tip Refill - Yellow Base (30)
S308	Intraoral Tip Refill - Dark Grey (50)



Blu-Mousse®

A Legacy VPS Bite Registration

- ▶ Stays stacked on the occlusal surface during placement.
- ▶ Optimal consistency allows patient to naturally close with zero resistance while capturing every detail.
- ▶ Hardens to a plaster-like final set with a durometer of 90 to ensure precise articulations.
- ▶ Set times: 30 Seconds | 60 Seconds | 2 Minutes

Additional Clinical Applications

- ▶ Scannable using the latest technologies.
- ▶ Used as a matrix for fabricating direct and indirect provisionals together with Mach Die Silicone by Parkell.



50 ML SPLIT-CARTRIDGE KITS:

Include (2) 50 ml cartridges (100 ml ttl.) and mixing tips.

SKU	DESCRIPTION	SET TIME	SCENT
S457S	Blu-Mousse Super-Fast	30 Sec.	Unscented
S448S	Blu-Mousse Super-Fast	30 Sec.	Peppermint
S440S	Blu-Mousse 60	60 Sec.	Unscented
S438S	Blu-Mousse Classic	2 Min.	Unscented
S437S	Blu-Mousse Classic	2 Min.	Peppermint

Accessories & Refills

S445	50 ml Split-Cartridge Mixing Tip Refill - Pink Base (30)
------	--

FROM IMPRESSION TO WHITENING TRAY

IN MINUTES! LEARN HOW TO DO IT.

Step 1 - Take alginate impression & inject Mach Die silicone.



Step 2 - Apply Blu-Mousse to achieve a stable base.



Step 3 - Separate the model from the impression.



Step 4 - Create a whitening tray ready for patient use.



Scan the QR code to watch the video tutorials