


AccuFilm® FASTCHECK

Pre-Mounted Articulating Film




Articulation. Simplified.

AccuFilm®, the world-renowned brand recognized for accurately marking occlusion for decades introduces its simplest and thinnest film to date. At only 8 microns, AccuFilm® FastCheck incorporates the accuracy dental professionals routinely rely on with the added convenience of a built-in frame.

- ✓ Built-in frame facilitates easy positioning over the occlusal surfaces of posterior teeth without the need for forceps
- ✓ Incredibly thin film at only 8 microns consistently produces sharp, easy-to-interpret markings even on hard-to-test highly polished metal and ceramic surfaces
- ✓ Easily identifies the smallest occlusal discrepancies with precise marking on the points of contact
- ✓ Double-sided color coating for simultaneous marking of both arches
- ✓ Equally effective on both wet and dry metal, porcelain, and zirconia surfaces



Ultra-thin! [_____]
Now only **8 Microns**



Item #	Product
S054	AccuFilm® FastCheck Double-Sided (Blue/Blue) Articulating Film; (100) Strips

See other side for our complete line of AccuFilm Articulating products →






Parkell Europe AB, Finvids väg 8, 194 27 Upplands Väsby, Sweden.

For more information, visit parkell.eu or contact infoeurope@parkell.com, Tel.+46 708 593 481.

parkell®




AccuFilm® Pre-Cut Film Canisters

Item #	Product
S014	 AccuFilm® Single-Sided Marking Film—Black; (1) canister of 280 pre-cut strips
S015	 AccuFilm® Single-Sided Marking Film—Red; (1) canister of 280 pre-cut strips
S017	 AccuFilm® Double-Sided Marking Film—Red/Black; (1) canister of 280 pre-cut strips
S018	 AccuFilm® Double-Sided Marking Film—Black/Black; (1) canister of 280 pre-cut strips
S021	 AccuFilm® Double-Sided Marking Film—Red/Red; (1) canister of 280 pre-cut strips



Single and double-sided film thickness is 19 microns.



AccuFilm® Booklet

Item #	Product
S053	 AccuFilm® Double-Sided Marking Film Booklet—Red/Black; (5) booklets (50 sheets/booklet, 250 sheets ttl.)

Double-sided film thickness is 19 microns.



Parkell Europe AB, Finvids väg 8, 194 27 Upplands Väsby, Sweden.

For more information, visit parkell.eu or contact infoeurope@parkell.com, Tel.+46 708 593 481.